

The Certain[®] Implant System Ordering Information

When It Clicks, It's Seated. You're Certain.



Tapered Implants

P = Restorative Platform
D = Diameter



OSSEOTITE® Tapered Certain PREVAIL

Commercially Pure Titanium

Length	4/3mm(P)	5/4mm(P)	6/5mm(P)
8.5mm	XIITP4385	XIITP5485	XIITP6585
10mm	XIITP4310	XIITP5410	XIITP6510
11.5mm	XIITP4311	XIITP5411	XIITP6511
13mm	XIITP4313	XIITP5413	XIITP6513
15mm	XIITP4315	XIITP5415	XIITP6515






3i T3 With DCD® Platform Switched

Commercially Pure Titanium

Length	4/3mm(P)	5/4mm(P)	6/5mm(P)
8.5mm	BNPT4385	BNPT5485	BNPT6585
10mm	BNPT4310	BNPT5410	BNPT6510
11.5mm	BNPT4311	BNPT5411	BNPT6511
13mm	BNPT4313	BNPT5413	BNPT6513
15mm	BNPT4315	BNPT5415	BNPT6515

Implant Cover Screws

	3.25mm(D)	4mm(D)	5mm(D)	6mm(D)
 Flat*	IMCSF34	ICSF41	ICSF50	ICSF60
 Standard	IMMCS1	ICS375	ICS500	ICS600
 Headless	N/A	ICS275	ICS275	ICS275

*All Implants include the corresponding Flat Cover Screw.

Tapered Implants

(Cont'd)

P = Restorative Platform
D = Diameter



Full OSSEOTITE Tapered Certain

Commercially Pure Titanium

Length	3.25mm(D)	4mm(D)	5mm(D)	6mm(D)
8.5mm	XIFNT3285	XIFNT485	XIFNT585	XIFNT685
10mm	XIFNT3210	XIFNT410	XIFNT510	XIFNT610
11.5mm	XIFNT3211	XIFNT411	XIFNT511	XIFNT611
13mm	XIFNT3213	XIFNT413	XIFNT513	XIFNT613
15mm	XIFNT3215	XIFNT415	XIFNT515	XIFNT615



OSSEOTITE® Tapered Certain

Titanium Alloy

Length	3.25mm(D)	4mm(D)	5mm(D)	6mm(D)
8.5mm	INT3285	INT485	INT585	INT685
10mm	INT3210	INT410	INT510	INT610
11.5mm	INT3211	INT411	INT511	INT611
13mm	INT3213	INT413	INT513	INT613
15mm	INT3215	INT415	INT515	INT615






3i T3 With DCD® Non-Platform Switched

Commercially Pure Titanium

Length	3.25mm(D)	4mm(D)	5mm(D)	6mm(D)
8.5mm	BNST3285	BNST485	BNST585	BNST685
10mm	BNST3210	BNST410	BNST510	BNST610
11.5mm	BNST3211	BNST411	BNST511	BNST611
13mm	BNST3213	BNST413	BNST513	BNST613
15mm	BNST3215	BNST415	BNST515	BNST615

Implant Cover Screws


	3.25mm(D)	4mm(D)	5mm(D)	6mm(D)
 Flat*	IMCSF34	ICSF41	ICSF50	ICSF60
 Standard	IMMCS1	ICS375	ICS500	ICS600
 Headless	N/A	ICS275	ICS275	ICS275

*All Implants include the corresponding Flat Cover Screw.


Parallel Walled Implants

(Cont'd)

P = Restorative Platform
D = Diameter




 **OSSEOTITE² Parallel Walled Certain**
 Commercially Pure Titanium

Length	3.25mm(D)	4mm(D)	5mm(D)	6mm(D)
8.5mm	XIFOSM385	XIFOSS485	XIFOSS585	XIFOSS685
10mm	XIFOSM310	XIFOSS410	XIFOSS510	XIFOSS610
11.5mm	XIFOSM311	XIFOSS411	XIFOSS511	XIFOSS611
13mm	XIFOSM313	XIFOSS413	XIFOSS513	XIFOSS613
15mm	XIFOSM315	XIFOSS415	XIFOSS515	XIFOSS615

 **OSSEOTITE² Parallel Walled Certain**
 Titanium Alloy

Length	3.25mm(D)	4mm(D)	5mm(D)	6mm(D)
8.5mm	IOSM385	IOSS485	IOSS585	IOSS685
10mm	IOSM310	IOSS410	IOSS510	IOSS610
11.5mm	IOSM311	IOSS411	IOSS511	IOSS611
13mm	IOSM313	IOSS413	IOSS513	IOSS613
15mm	IOSM315	IOSS415	IOSS515	IOSS615
18mm	IOSM318	IOSS418	N/A	N/A
20mm	N/A	IOSS420	N/A	N/A

Implant Cover Screws





	3.25mm(D)	4mm(D)	5mm(D)	6mm(D)
 Flat*	IMCSF34	ICSF41	ICSF50	ICSF60
 Standard	IMMCS1	ICS375	ICS500	ICS600
 Headless	N/A	ICS275	ICS275	ICS275


*All Implants include the corresponding Flat Cover Screw.

Healing Abutments (Cont'd)

D = Diameter
P = Emergence Profile
H = Collar Height

EP® One-Piece Healing Abutments


Dimensions	Item #	Dimensions	Item #
3.4mm(D) 	3.8mm(P) X 2mm(H)	5mm(D) 	5.6mm(P) X 2mm(H)
	3.8mm(P) X 3mm(H)		5.6mm(P) X 3mm(H)
	3.8mm(P) X 4mm(H)		5.6mm(P) X 4mm(H)
	3.8mm(P) X 6mm(H)		5.6mm(P) X 6mm(H)
			5.6mm(P) X 8mm(H)
	5mm(P) X 2mm(H)		6mm(P) X 2mm(H)
	5mm(P) X 3mm(H)		6mm(P) X 3mm(H)
	5mm(P) X 4mm(H)		6mm(P) X 4mm(H)
	5mm(P) X 6mm(H)		6mm(P) X 6mm(H)
			6mm(P) X 8mm(H)
			7.5mm(P) X 2mm(H)
			7.5mm(P) X 3mm(H)
	7.5mm(P) X 4mm(H)		
	7.5mm(P) X 6mm(H)		
	7.5mm(P) X 8mm(H)		
4.1mm(D) 	5mm(P) X 2mm(H)	6mm(D) 	6.8mm(P) X 2mm(H)
	5mm(P) X 3mm(H)		6.8mm(P) X 3mm(H)
	5mm(P) X 4mm(H)		6.8mm(P) X 4mm(H)
	5mm(P) X 6mm(H)		6.8mm(P) X 6mm(H)
	5mm(P) X 8mm(H)		6.8mm(P) X 8mm(H)
	6mm(P) X 3mm(H)		7.5mm(P) X 2mm(H)
	6mm(P) X 4mm(H)		7.5mm(P) X 3mm(H)
	6mm(P) X 6mm(H)		7.5mm(P) X 4mm(H)
	6mm(P) X 8mm(H)		7.5mm(P) X 6mm(H)
			7.5mm(P) X 8mm(H)
	7.5mm(P) X 3mm(H)		
	7.5mm(P) X 4mm(H)		
	7.5mm(P) X 6mm(H)		
	7.5mm(P) X 8mm(H)		




Impression Copings, Laboratory Components & Provisional Restorations

Impression Copings

Pick-Up Impression Coping


	Description	Item #
	3.4mm(D) 3.4mm(P)	IMIC30
	3.8mm(P)	IMIC33
	5mm(P)	IMIC35
4.1mm(D)	4.1mm(P)	IIIC41
	4.1mm(P) (Non-Hexed)	IIIC42
	5mm(P)	IIIC12
	6mm(P)	IIIC60
	7.5mm(P)	IIIC75
5mm(D)	5mm(P)	IWIP50
	5.6mm(P)	IWIP55
	6mm(P)	IWIP56
	7.5mm(P)	IWIP57
6mm(D)	6mm(P)	IWIP60
	6.8mm(P)	IWIP66
	7.5mm(P)	IWIP67

Twist Lock™ Transfer Impression Coping


	Description	Item #
	3.4mm(D) 3.8mm(P)	IMIT33
	5mm(P)	IMIT35
4.1mm(D)	4.1mm(P)	IIIC44
	5mm(P)	IIIC45
	6mm(P)	IIIC46
	7.5mm(P)	IIIC47
5mm(D)	5mm(P)	IWIT50
	5.6mm(P)	IWIT55
	6mm(P)	IWIT56
	7.5mm(P)	IWIT57
6mm(D)	6mm(P)	IWIT60
	6.8mm(P)	IWIT66
	7.5mm(P)	IWIT67

Laboratory Components

Implant Laboratory Analog

Description	Item #
	
3.4mm(D)	IMMILA
4.1mm(D)	IILA20
5mm(D)	IILAW5
6mm(D)	IILAW6

Titanium Temporary Cylinder*

Description	Item #	Item #
		
	Hexed	Non-Hexed
3.4mm(D)	IMMTCS1	N/A
4.1mm(D)	IITCS41	IITCS42
5mm(D)	IWTCS51	IWTCS52
6mm(D)	IWTCS61	IWTCS62


D = Diameter
P = Emergence Profile
H = Collar Height

Cement-Retained Restorations

(Cont'd)

D = Diameter
P = Emergence Profile
H = Collar Height

GingiHue® Post*





		Description	Item #	Item #
			Straight	Pre-Angled (15°)
	3.4mm(D)	3.8mm(P) X 2mm(H)	IMAP32G	IMPAP32G
		3.8mm(P) X 4mm(H)	IMAP34G	IMPAP34G
		5mm(P) X 2mm(H)	IMAP352G	IMPAP352G
		5mm(P) X 4mm(H)	IMAP354G	IMPAP354G
4.1mm(D)	5mm(P) X 2mm(H)	IAPP452G	IPAP452G	
	5mm(P) X 4mm(H)	IAPP454G	IPAP454G	
	6mm(P) X 2mm(H)	IAPP462G	IPAP462G	
	6mm(P) X 4mm(H)	IAPP464G	IPAP464G	
	7.5mm(P) X 2mm(H)	IAPP472G	IPAP472G	
	7.5mm(P) X 4mm(H)	IAPP474G	IPAP474G	
5mm(D)	5mm(P) X 2mm(H)	IWPP552G	IPAP552G	
	5mm(P) X 4mm(H)	IWPP554G	IPAP554G	
	6mm(P) X 2mm(H)	IWPP562G	IPAP562G	
	6mm(P) X 4mm(H)	IWPP564G	IPAP564G	
	7.5mm(P) X 2mm(H)	IWPP572G	IPAP572G	
	7.5mm(P) X 4mm(H)	IWPP574G	IPAP574G	
6mm(D)	6mm(P) X 2mm(H)	IWPP662G	IPAP662G	
	6mm(P) X 4mm(H)	IWPP664G	IPAP664G	
	7.5mm(P) X 2mm(H)	IWPP672G	IPAP672G	
	7.5mm(P) X 4mm(H)	IWPP674G	IPAP674G	

*GingiHue Posts do not include the screw. The final Gold-Tite® Screw (IUNIHG) must be ordered separately.

Cement-Retained Restorations (Cont'd)

D = Diameter

UCLA Abutment* (For Use With Screw-Retained Restorations Unless Otherwise Noted)

	Description		Item #
	3.4mm(D)	Gold Hexed	IMUCG1C* (Cement- Or Screw-Retained)
		Gold Non-Hexed (With Gold-Tite® Large Diameter Screw)	IMUCG2C
		Gold Non-Hexed (With Ti Large Diameter Screw)	IMUCG2T
		Castable Non-Hexed (With Gold-Tite Large Diameter Screw)	IMUCC2C
		Castable Non-Hexed (With Ti Large Diameter Screw)	IMUCC2T
	4.1mm(D)	Gold Hexed	IGUCA1C* (Cement- Or Screw-Retained)
		Gold Non-Hexed (With Gold-Tite Large Diameter Screw)	IGUCA2C
		Gold Non-Hexed (With Ti Large Diameter Screw)	IGUCA2T
		Castable Non-Hexed (With Gold-Tite Large Diameter Screw)	IUNAB2C
		Castable Non-Hexed (With Ti Large Diameter Screw)	IUNAB2T
	5mm(D)	Gold Hexed	IWGA51C* (Cement- Or Screw-Retained)
		Gold Non-Hexed (With Gold-Tite Large Diameter Screw)	IWGA52C
		Gold Non-Hexed (With Ti Large Diameter Screw)	IWGA52T
		Castable Non-Hexed (With Gold-Tite Large Diameter Screw)	IWPC52C
		Castable Non-Hexed (With Ti Large Diameter Screw)	IWPC52T
	6mm(D)	Gold Hexed	IWGA61C* (Cement- Or Screw-Retained)
		Gold Non-Hexed (With Gold-Tite Large Diameter Screw)	IWGA62C
		Gold Non-Hexed (With Ti Large Diameter Screw)	IWGA62T
		Castable Non-Hexed (With Gold-Tite Large Diameter Screw)	IWPC62C
		Castable Non-Hexed (With Ti Large Diameter Screw)	IWPC62T

*UCLA Hexed Abutments do not include the screw. The final Gold-Tite Screw (IUNIHG) or the Titanium Screw (IUNIHT) must be ordered separately.

Screw-Retained Restorations




Low Profile Abutment*

	Collar Height	One Piece	Two Piece	17°	30°
Seating Surface 3.4mm(D)	1mm	ILPC341U	ILPC341	-	-
	2mm	ILPC342U	ILPC342	ILPAC3217	-
	3mm	ILPC343U	ILPC343	-	ILPAC3330
	4mm	ILPC344U	ILPC344	ILPAC3417	-
	5mm	-	-	-	ILPAC3530
4.1mm(D)	1mm	ILPC441U	ILPC441	-	-
	2mm	ILPC442U	ILPC442	ILPAC4217	-
	3mm	ILPC443U	ILPC443	-	ILPAC4330
	4mm	ILPC444U	ILPC444	ILPAC4417	-
	5mm	-	-	-	ILPAC4530
5mm(D)	1mm	ILPC541U	ILPC541	-	-
	2mm	ILPC542U	ILPC542	ILPAC5217	-
	3mm	ILPC543U	ILPC543	-	ILPAC5330
	4mm	ILPC544U	ILPC544	ILPAC5417	-
	5mm	-	-	-	ILPAC5530



*Requires the use of an Abutment Driver (PAD00, PAD24) or Driver Tip (RASA3).

Low Profile Screws


Description	Item #
 Gold-Tite® Retaining Screw	LPCGSH*
Titanium Retaining Screw	LPCTSH*
Waxing Screw	LPCWS

*The recommended torque for the Retaining Screws is 10Ncm.



Low Profile Laboratory Tools

Description	Item #
 Lapping Tool	LPCAMI
 Polishing Protector	LPCPP


Low Profile Impression Copings & Analog

Description	Item #
 Pick Up Impression Coping (Hexed)	LPCPIC1
Pick Up Impression Coping (Non-Hexed)	LPCPIC2
Twist Lock™ Impression Coping (Hexed)	LPCTIC1
Twist Lock Impression Coping (Non-Hexed)	LPCTIC2
QuickBridge® Impression Coping	LPCRIC
Laboratory Analog	LPCLA


Low Profile Provisional Restoration Components

Description	Item #
 Temporary Cylinder (Hexed)	LPCTC1
Temporary Cylinder (Non-Hexed)	LPCTC2
PreFormance® Temporary Cylinder (Hexed)	LPCPTC1
PreFormance Temporary Cylinder (Non-Hexed)	LPCPTC2
 Healing Cap	LPCHC
QuickBridge	LPCQB
QuickBridge Replacement Cap	LPCQBCAP
Distal Extension	LPCDE

Low Profile Final Restoration Components

Description	Item #
 Gold Cylinder (Hexed)	LPCGC1
Gold Cylinder (Non-Hexed)	LPCGC2
Castable Cylinder (Hexed)	LPCCC1
Castable Cylinder (Non-Hexed)	LPCCC2
Laser Welded Cylinder	LPCTUCA


Large Hexed Driver & Driver Tip

Description	Item #
 Large Hexed Driver	PHD02N, PHD03N*
Large Hexed Driver Tip	RASH3N, RASH8N*
Posterior Abutment Driver 17mm(L)	PAD00
Standard Abutment Driver 24mm(L)	PAD24

*Made narrower to fit with external hexed angled Low Profile Abutments. The Driver and Driver Tips can be identified by a laser marked dot after the item number as pictured below.







Low Profile Bone Profilers


Description	Item #
 Bone Profiler, 3.4mm(D)	LPCBP3
Bone Profiler, 4.1mm(D)	LPCBP4
Bone Profiler, 5mm(D)	LPCBP5

Overdenture Restorations















LOCATOR® Abutment*

	Description	Item #
Seating Surface  	1mm(H)	IMLOA001
	2mm(H)	IMLOA002
	3mm(H)	IMLOA003
	4mm(H)	IMLOA004
	5mm(H)	IMLOA005
	6mm(H)	IMLOA006
 	1mm(H)	ILOA001
	2mm(H)	ILOA002
	3mm(H)	ILOA003
	4mm(H)	ILOA004
	5mm(H)	ILOA005
	6mm(H)	ILOA006

*Included with each abutment:








LOCATOR Abutment Components

Description	Item #
 Impression Coping (one time use)	LAIC1
 Laboratory Analog	LALA1
 Core Tool/Abutment Driver	LCTDR1
 Driver Tip 24mm/30mm	LOADT4/LOADT9
 Replacement Housing Kit	LORHK
 Extended Range Males (4-pack) - Green	LAERM
 Extra Light Retention Males (4-pack) - Blue	LAELM
 Light Retention Replacement Males (4-pack) - Pink	LLRMS
 Replacement Males (4-pack) - White	LARMS
 Extra Light Retention Extended Range Males - Red	LELARM
 Replacement Housing	LOAH
 1.25mm .050in. Driver Tip 24mm(L) - used with LCTDR1	RASH4
 1.25mm .050in. Driver Tip 30mm(L) - used with LCTDR1	RASH9
 Male Extraction Tool	LOA8397*

*Part of the Core Tool/Abutment Driver LCTDR1.





Screws, Drivers, Driver Tips & Misc. Items

Abutment Screws (UCLA, GingiHue®, Provide®)

Description	Item #
 Gold-Tite® Hexed Screw	IUNIHG*
 Titanium Hexed Screw	IUNIHT*
 Gold-Tite Hexed Large Diameter Screw	ILRGHG*
 Titanium Hexed Large Diameter Screw	ILRGHT*
 ZiReal® Gold-Tite Hexed Screw	IZSHG*




*The recommended torque for the Abutment Screws is 20Ncm.

Try-In & Impression Coping Screws

Description	Item #
 Hexed Try-In Screw (5-pack)	IUNITS*
 Knurled Implant Waxing Screw/Guide Pin	IWSU30
 Twist Lock™ Coping Screw	IWITSC
 ZiReal Titanium Try-in Screw	IZIRTS

*Only for GingiHue Posts, Hexed UCLA Abutments and Hexed Titanium Temporary Cylinders.

Drivers, Driver Tips & Miscellaneous Items

Description	Item #
 Narrow Posterior Large Hexed Driver 17mm	PHD02N
 Narrow Standard Large Hexed Driver 24mm	PHD03N
 Posterior Abutment Driver 17mm	PAD00
 Posterior Abutment Driver 24mm	PAD24
 Manual Interchangeable Driver Tip Handle	MIDTH
 Twist Lock Coping Driver	ICD00
 Laboratory Abutment Holder (all diameters)	ILTAH57
 Certain® QuickSeat® Activator Tool	IQSA01
 Screw Removal Kit	ISRT10
 Narrow Right Angle Hexed Driver Tip 24mm	RASH3N
 Narrow Right Angle Hexed Driver Tip 30mm	RASH8N
 Right Angle Abutment Driver Tip - Steel	RASA3

D = Diameter
H = Collar Height

Tapered Implants Surgical Kits



P = Restorative Platform
D = Diameter
L = Length

QNTSK Contents

QNTSK40 Premium Kit

QNTSK40U Upgrade Kit

QNTSK20 Starter Kit

QNTSK40 Premium Kit

QNTSK40U Upgrade Kit

QNTSK20 Starter Kit

Description	Item #	Description	Item #
Tapered Implant Depth Indicator 3.25mm(D) X 8.5mm(L) 3 Ea.	NTDI3285	Quad Shaping Drill 6mm(D) X 13mm(L)	QSD613
Tapered Implant Depth Indicator 3.25mm(D) X 10mm(L)	NTDI3210	Quad Shaping Drill 6mm(D) X 15mm(L)	QSD615
Tapered Implant Depth Indicator 3.25mm(D) X 11.5mm(L)	NTDI3211	Narrow Posterior Hexed Driver - 17mm(L)	PHD00N
Tapered Implant Depth Indicator 3.25mm(D) X 13mm(L)	NTDI3213	Posterior Large Hexed Driver - 17mm(L)	PHD02N
Tapered Implant Depth Indicator 3.25mm(D) X 15mm(L)	NTDI3215	Narrow Right Angle Large Hexed Driver Tip - 24mm(L)	RASH3N
Tapered Implant Depth Indicator 4mm(D) X 8.5mm(L) 3 Ea.	NTDI485	Certain Universal Ratchet Ext. Driver 4.1, 5, 6mm(D) - Short	IRE100U
Tapered Implant Depth Indicator 4mm(D) X 10mm(L)	NTDI410	Certain Universal Ratchet Ext. Driver 4.1, 5, 6mm(D) - Long	IRE200U
Tapered Implant Depth Indicator 4mm(D) X 11.5mm(L)	NTDI411	Direction Indicator 2mm(D) X 10mm(L) 3 Ea.	DI100
Tapered Implant Depth Indicator 4mm(D) X 13mm(L)	NTDI413	ACT Drill Extension	ACTDE
Tapered Implant Depth Indicator 4mm(D) X 15mm(L)	NTDI415	Handpiece Connector	MDR10
Tapered Implant Depth Indicator 5mm(D) X 8.5mm(L) 3 Ea.	NTDI585	Implant Mount 3.4mm(D) X 15mm(L)	MMC15
Tapered Implant Depth Indicator 5mm(D) X 10mm(L)	NTDI510	Certain Implant Driver Tip 4, 5, 6mm(D) - Long	IIPDTL
Tapered Implant Depth Indicator 5mm(D) X 11.5mm(L)	NTDI511	Certain Implant Driver Tip 4, 5, 6mm(D) - Short	IIPDTS
Tapered Implant Depth Indicator 5mm(D) X 13mm(L)	NTDI513	Certain Implant Driver Tip 3.4mm(D) - Short	IMPDTS
Tapered Implant Depth Indicator 5mm(D) X 15mm(L)	NTDI515	Certain Implant Driver Tip 3.4mm(D) - Long	IMPDTL
Tapered Implant Depth Indicator 6mm(D) X 8.5mm(L) 3 Ea.	NTDI685	Tapered Implant Dense Bone Tap 3.25mm(D) X 8.5mm(L)	NTAP3285
Tapered Implant Depth Indicator 6mm(D) X 10mm(L)	NTDI610	Tapered Implant Dense Bone Tap 3.25mm(D) X 10mm(L)	NTAP3210
Tapered Implant Depth Indicator 6mm(D) X 11.5mm(L)	NTDI611	Tapered Implant Dense Bone Tap 3.25mm(D) X 11.5mm(L)	NTAP3211
Tapered Implant Depth Indicator 6mm(D) X 13mm(L)	NTDI613	Tapered Implant Dense Bone Tap 3.25mm(D) X 13mm(L)	NTAP3213
Tapered Implant Depth Indicator 6mm(D) X 15mm(L)	NTDI615	Tapered Implant Dense Bone Tap 3.25mm(D) X 15mm(L)	NTAP3215
ACT® Pointed Starter Drill	ACTPSD	Tapered Implant Dense Bone Tap 4mm(D) X 8.5mm(L)	NTAP485
ACT Reusable Twist Drill 2mm(D) X 10mm(L)	ACT2010	Tapered Implant Dense Bone Tap 4mm(D) X 10mm(L)	NTAP410
ACT Reusable Twist Drill 2mm(D) X 15mm(L)	ACT2015	Tapered Implant Dense Bone Tap 4mm(D) X 11.5mm(L)	NTAP411
ACT Reusable Twist Drill 2mm(D) X 20mm(L)	ACT2020	Tapered Implant Dense Bone Tap 4mm(D) X 13mm(L)	NTAP413
Certain® Countersink Drill, PREVAIL® 3/4/3mm(P)	ICD100	Tapered Implant Dense Bone Tap 4mm(D) X 15mm(L)	NTAP415
Quad Shaping Drill 3.25mm(D) X 8.5mm(L)	QSD3285	Tapered Implant Dense Bone Tap 5mm(D) X 8.5mm(L)	NTAP585
Quad Shaping Drill 3.25mm(D) X 10mm(L)	QSD3210	Tapered Implant Dense Bone Tap 5mm(D) X 10mm(L)	NTAP510
Quad Shaping Drill 3.25mm(D) X 11.5mm(L)	QSD3211	Tapered Implant Dense Bone Tap 5mm(D) X 11.5mm(L)	NTAP511
Quad Shaping Drill 3.25mm(D) X 13mm(L)	QSD3213	Tapered Implant Dense Bone Tap 5mm(D) X 13mm(L)	NTAP513
Quad Shaping Drill 3.25mm(D) X 15mm(L)	QSD3215	Tapered Implant Dense Bone Tap 5mm(D) X 15mm(L)	NTAP515
Quad Shaping Drill 4mm(D) X 8.5mm(L)	QSD485	Tapered Implant Dense Bone Tap 6mm(D) X 8.5mm(L)	NTAP685
Quad Shaping Drill 4mm(D) X 10mm(L)	QSD410	Tapered Implant Dense Bone Tap 6mm(D) X 10mm(L)	NTAP610
Quad Shaping Drill 4mm(D) X 11.5mm(L)	QSD411	Tapered Implant Dense Bone Tap 6mm(D) X 11.5mm(L)	NTAP611
Quad Shaping Drill 4mm(D) X 13mm(L)	QSD413	Tapered Implant Dense Bone Tap 6mm(D) X 13mm(L)	NTAP613
Quad Shaping Drill 4mm(D) X 15mm(L)	QSD415	Tapered Implant Dense Bone Tap 6mm(D) X 15mm(L)	NTAP615
Quad Shaping Drill 5mm(D) X 8.5mm(L)	QSD585	Open End Wrench	CW100
Quad Shaping Drill 5mm(D) X 10mm(L)	QSD510	Certain Driver Replacement O-Rings (10-pack)	IRORDR
Quad Shaping Drill 5mm(D) X 11.5mm(L)	QSD511	Tapered Implant Dense Bone Tap Kit Tray	NTAPKT
Quad Shaping Drill 5mm(D) X 13mm(L)	QSD513	Ratchet Extension Short 6mm(L)	RE100
Quad Shaping Drill 5mm(D) X 15mm(L)	QSD515	Ratchet Extension Long 15mm(L)	RE200
Quad Shaping Drill 6mm(D) X 8.5mm(L)	QSD685	Ratchet Wrench	WR150
Quad Shaping Drill 6mm(D) X 10mm(L)	QSD610	Base, Molded Insert & Smoked Lid	QNTSKT
Quad Shaping Drill 6mm(D) X 11.5mm(L)	QSD611	Implant Depth Gauge 20mm(L)	DPO20

Parallel Walled Implants Surgical Kits



D = Diameter
L = Length

PSKT Contents

PSKT40 Premium Kit

PSKT35 Certain® Premium Kit

PSKT01 Starter Kit

PSKT40 Premium Kit

PSKT35 Certain Premium Kit

PSKT01 Starter Kit

Description	Item #				Description	Item #			
Plastic Surgical Tray Organizer	PTT300I	•	•	•	Countersink Depth Indicator for CD5600	CDI5600	•	•	•
Plastic Tray Organizer for PTT300I	PTT100				Certain Countersink Drill PREVAIL 3/4/3	ICD100	•	•	•
ACT® Reusable Twist Drill 2mm(D) X 10mm(L)	ACT2010				Countersink Depth Indicator for ICD100	ICDI100	•	•	•
ACT Reusable Twist Drill 2mm(D) X 15mm(L)	ACT2015	•	•	•	Countersink Drill 4.1mm	CD100	•		•
ACT Reusable Twist Drill 2mm(D) X 20mm(L)	ACT2020				Countersink/Pilot Drill 6mm(D)	CD600		•	•
ACT Reusable Twist Drill 2.75mm(D) X 10mm(L)	ACT2710				Pilot Drill 2-3mm(D)	PD100	•	•	•
ACT Reusable Twist Drill 2.75mm(D) X 15mm(L)	ACT2715	•	•	•	Implant Dense Bone Tap 3.25mm(D) X 18mm(L)	MTAP2		•	•
ACT Reusable Twist Drill 2.75mm(D) X 20mm(L)	ACT2720				Implant Dense Bone Tap 4mm(D) X 13mm(L)	TAP413		•	•
ACT Reusable Twist Drill 3mm(D) X 10mm(L)	ACT3010				Implant Dense Bone Tap 5mm(D) X 13mm(L)	TAP53S		•	•
ACT Reusable Twist Drill 3mm(D) X 15mm(L)	ACT3015	•	•	•	Implant Dense Bone Tap 6mm(D) X 13mm(L)	TAP63S			•
ACT Reusable Twist Drill 3mm(D) X 20mm(L)	ACT3020				OSSEOTITE® 2 Dense Bone Tap 5mm(D) X 7-13mm(L)	XTAP53S		•	•
ACT Reusable Twist Drill 3.25mm(D) X 10mm(L)	ACT3210				OSSEOTITE 2 Dense Bone Tap 6mm(D) X 7-13mm(L)	XTAP63S			•
ACT Reusable Twist Drill 3.25mm(D) X 15mm(L)	ACT3215	•	•	•	Certain Implant Driver Tip (Short) 3.4mm(D)	IMPPTS		•	•
ACT Reusable Twist Drill 3.25mm(D) X 20mm(L)	ACT3220				Certain Implant Driver Tip (Long) 3.4mm(D)	IMPDTL	•	•	•
ACT Reusable Twist Drill 3.85mm(D) X 10mm(L)	ACT3810	•	•	•	Certain Implant Placement Driver Tip (Short) 4.1, 5, 6mm(D)	IIPPTS	•	•	•
ACT Reusable Twist Drill 3.85mm(D) X 15mm(L)	ACT3815	•	•	•	Certain Implant Placement Driver Tip (Long) 4.1, 5, 6mm(D)	IIPDTL	•	•	•
ACT Reusable Twist Drill 3.85mm(D) X 20mm(L)	ACT3820	•	•	•	Certain Universal Ratchet Extension Driver (Short) 6mm(L)	IRE100U	•	•	•
ACT Reusable Twist Drill 4.25mm(D) X 8.5mm(L)	ACT4285				Certain Universal Ratchet Extension Driver (Long) 5mm(L)	IRE200U	•	•	•
ACT Reusable Twist Drill 4.25mm(D) X 13mm(L)	ACT4213	•	•	•	Narrow Small Hex Hand Driver 17mm(L)	PHD00N	•		•
ACT Reusable Twist Drill 4.25mm(D) X 18mm(L)	ACT4218				Narrow Large Hex Hand Driver 17mm(L)	PHD02N	•	•	•
ACT Reusable Twist Drill 4.85mm(D) X 8.5mm(L)	ACT4885				Narrow Large Hex Driver Tip 24mm(L)	RASH3N	•	•	•
ACT Reusable Twist Drill 4.85mm(D) X 13mm(L)	ACT4813	•	•	•	Driver Tip Replacement O-Rings (10-pack)	IRORDR		•	•
ACT Reusable Twist Drill 4.85mm(D) X 18mm(L)	ACT4818				Handpiece Connector	MDR10	•	•	•
ACT Reusable Twist Drill 5.25mm(D) X 8.5mm(L)	ACT5285	•	•	•	Open End Wrench	CW100	•		•
ACT Reusable Twist Drill 5.25mm(D) X 13mm(L)	ACT5213				Ratchet Wrench	WR150	•	•	•
ACT Reusable Twist Drill 5.25mm(D) X 18mm(L)	ACT5218				Ratchet Extension (Short) 6mm(L)	RE100	•		•
ACT Drill Extension	ACTDE	•	•	•	Ratchet Extension (Long) 15mm(L)	RE200	•	•	•
Countersink Drill OSSEOTITE XP® 4/5 PREVAIL® 4/5/4	CD4500				Implant Depth Gauge 20mm(L)	DP020	•	•	•
Countersink Depth Indicator for CD4500	CDI4500	•	•	•	Implant Mount 3.4mm(D) X 15mm(L)	MMC15	•		•
Countersink/Pilot Drill 5mm(D)	CD500	•	•	•	Direction Indicator 2mm(D) X 10mm(L) 3ea.	DI100	•	•	•
Countersink Drill OSSEOTITE XP 5/6 PREVAIL 5/6/5	CD5600				Stainless Steel Round Drill 2mm(D)	RD100	•	•	•

BIOMET **3i** offers a family of synergistic system solutions to help you and your patients achieve beautiful, confident and healthy smiles.



■ Smile Beautifully ■ Smile Confidently ■ Smile Healthy

For more information regarding Smile Therapy™ System Solutions, please contact your local BIOMET **3i** Sales Representative.

BIOMET 3i™

Check Out The BIOMET **3i** Online Store at www.shopbiomet3i.com
(Available for Spain & US Customers Only)



Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Toll Free: 1-800-342-5454
Phone: +1-561-776-6700
Fax: +1-561-776-1272
www.biomet3i.com

EC REP BIOMET **3i** Dental Ibérica S.L.
EMEA Headquarters
WTC Almeda Park, Ed. 4, Planta 2
C/Tirso de Molina, 40
08940, Cornellà de Llobregat (Barcelona) Spain
Phone: +34-93-470-55-00
Fax: +34-93-371-78-49
www.biomet3i.com



All trademarks herein are the property of BIOMET **3i** LLC unless otherwise indicated. LOCATOR is a registered trademark of Zest Anchors, Inc. ©2015 BIOMET **3i** LLC.

This material is intended for clinicians only and is NOT intended for patient distribution. This material is not to be redistributed, duplicated or disclosed without the express written consent of BIOMET **3i**. For additional product information, including indications, contraindications, warnings, precautions and potential adverse effects, please visit the BIOMET **3i** Website: www.ifu.biomet3i.com.



ART1048
REV H 11/15