

Gold-Tite[®] Screw

Preservation By Design[®]




Implant-Abutment Seal Strength with the BIOMET **3i** Gold-Tite® Screw

A stable, tight implant/abutment interface minimizes abutment micromotion and reduces potential microleakage.

- Use of a BIOMET **3i** Certain® Gold-Tite Screw increases implant/abutment clamping force by 113% versus a non-coated screw.^{1*}
- The patented** Gold-Tite surface lubrication allows the screw to rotate further, increasing the clamping force and maximizing abutment stability.

*Pre-clinical studies are not necessarily indicative of clinical results.



Screw		Driver			Torque
Certain® Gold-Tite Hexed Screw	IUNIHG	Hexed Hand Driver	PHD02N	17mm(L)	20Ncm
			PHD03N	24mm(L)	
		Hexed Driver Tip	RASH3N	25mm(L)	
			RASH8N	30mm(L)	

¹Suttin Z¹, Towse R¹. Effect of abutment screw design in implant system seal performance. Presented at the European Association for Osseointegration, 20th Annual Scientific Meeting; October 2012; Copenhagen, Denmark. http://biomet3i.com/Pdf/Posters/P-450_Effect_of_Screw_Design_on_Implant_Seal_pdf

[†]The authors conducted this research while employed at BIOMET **3i**.

**USPTO 6287116 and additional patents pending.

*For More Information, Please Contact Your Local BIOMET **3i** Sales Representative.*

www.biomet3i.com



Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
1-800-342-5454
Outside the U.S.: +1-561-776-6700
Fax: +1-561-776-1272
www.biomet3i.com



BIOMET **3i**
Dental Iberica S.L.
WTC Alameda Park, Ed. 1, Planta 1ª
Pl. de la Pau, s/n
08940, Cornellà de Llobregat
(Barcelona) Spain
Phone: +34-93-470-55-00
Fax: +34-93-371-78-49



Join Us



Follow Us



Watch Us



Download It

Certain, Gold-Tite and Preservation By Design are registered trademarks and Providing Solutions – One Patient At A Time and design are trademarks of BIOMET **3i** LLC. ©2014 BIOMET **3i** LLC.

All trademarks herein are the property of BIOMET **3i** LLC unless otherwise indicated. This material is intended for clinicians only and is NOT intended for patient distribution. This material is not to be redistributed, duplicated, or disclosed without the express written consent of BIOMET **3i**. For additional product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see the product package insert and the BIOMET **3i** Website. For instructions go to: ifu.biomet3i.com.



INST1249
REV A 04/14