


QNTSK Surgical Kits & Components Assembly

Simplicity By Design




The Tapered Implant Surgical Tray –

(NOTE: Below Graphic Illustrates Component




Simplicity In Site Preparation

(Location For QNTSK40 Reference)



Tapered Implant Dense Bone Tap Kit



25

NTAP32XX



26

NTAP4XX



27

NTAP5XX



28

NTAP6XX

NOTE: Tapping in dense bone is required. Exceeding an insertion torque of more than 90Ncm may deform or strip the driver tip or the implant's internal hex and may possibly delay the surgical procedure.

QNTSK Surgical Tray Contents

QNTSK40 Premium Kit

QNTSK40U Upgrade Kit

QNTSK20 Starter Kit

QNTSK40 Premium Kit

QNTSK40U Upgrade Kit

QNTSK20 Starter Kit

Description		Item #		
1	Tapered Implant Depth Indicator 3.25mm(D) X 8.5mm(L) 3 Ea.	NTDI3285	•	•
1	Tapered Implant Depth Indicator 3.25mm(D) X 10mm(L)	NTDI3210	•	•
1	Tapered Implant Depth Indicator 3.25mm(D) X 11.5mm(L)	NTDI3211	•	•
1	Tapered Implant Depth Indicator 3.25mm(D) X 13mm(L)	NTDI3213	•	•
1	Tapered Implant Depth Indicator 3.25mm(D) X 15mm(L)	NTDI3215	•	•
2	Tapered Implant Depth Indicator 4mm(D) X 8.5mm(L) 3 Ea.	NTDI485	•	•
2	Tapered Implant Depth Indicator 4mm(D) X 10mm(L)	NTDI410	•	•
2	Tapered Implant Depth Indicator 4mm(D) X 11.5mm(L)	NTDI411	•	•
2	Tapered Implant Depth Indicator 4mm(D) X 13mm(L)	NTDI413	•	•
2	Tapered Implant Depth Indicator 4mm(D) X 15mm(L)	NTDI415	•	•
3	Tapered Implant Depth Indicator 5mm(D) X 8.5mm(L) 3 Ea.	NTDI585	•	•
3	Tapered Implant Depth Indicator 5mm(D) X 10mm(L)	NTDI510	•	•
3	Tapered Implant Depth Indicator 5mm(D) X 11.5mm(L)	NTDI511	•	•
3	Tapered Implant Depth Indicator 5mm(D) X 13mm(L)	NTDI513	•	•
3	Tapered Implant Depth Indicator 5mm(D) X 15mm(L)	NTDI515	•	•
4	Tapered Implant Depth Indicator 6mm(D) X 8.5mm(L) 3 Ea.	NTDI685	•	•
4	Tapered Implant Depth Indicator 6mm(D) X 10mm(L)	NTDI610	•	•
4	Tapered Implant Depth Indicator 6mm(D) X 11.5mm(L)	NTDI611	•	•
4	Tapered Implant Depth Indicator 6mm(D) X 13mm(L)	NTDI613	•	•
4	Tapered Implant Depth Indicator 6mm(D) X 15mm(L)	NTDI615	•	•
5	ACT® Pointed Starter Drill	ACTPSD	•	•
6	ACT Reusable Twist Drill 2mm(D) X 10mm(L)	ACT2010	•	•
6	ACT Reusable Twist Drill 2mm(D) X 15mm(L)	ACT2015	•	•
6	ACT Reusable Twist Drill 2mm(D) X 20mm(L)	ACT2020	•	•
7	Certain® Countersink Drill, PREVAIL® 3/4/3mm(P)	ICD100	•	•
8	Quad Shaping Drill 3.25mm(D) X 8.5mm(L)	QSD3285	•	•
8	Quad Shaping Drill 3.25mm(D) X 10mm(L)	QSD3210	•	•
8	Quad Shaping Drill 3.25mm(D) X 11.5mm(L)	QSD3211	•	•
8	Quad Shaping Drill 3.25mm(D) X 13mm(L)	QSD3213	•	•
8	Quad Shaping Drill 3.25mm(D) X 15mm(L)	QSD3215	•	•
9	Quad Shaping Drill 4mm(D) X 8.5mm(L)	QSD485	•	•
9	Quad Shaping Drill 4mm(D) X 10mm(L)	QSD410	•	•
9	Quad Shaping Drill 4mm(D) X 11.5mm(L)	QSD411	•	•
9	Quad Shaping Drill 4mm(D) X 13mm(L)	QSD413	•	•
9	Quad Shaping Drill 4mm(D) X 15mm(L)	QSD415	•	•
10	Quad Shaping Drill 5mm(D) X 8.5mm(L)	QSD585	•	•
10	Quad Shaping Drill 5mm(D) X 10mm(L)	QSD510	•	•
10	Quad Shaping Drill 5mm(D) X 11.5mm(L)	QSD511	•	•
10	Quad Shaping Drill 5mm(D) X 13mm(L)	QSD513	•	•
10	Quad Shaping Drill 5mm(D) X 15mm(L)	QSD515	•	•
11	Quad Shaping Drill 6mm(D) X 8.5mm(L)	QSD685	•	•
11	Quad Shaping Drill 6mm(D) X 10mm(L)	QSD610	•	•
11	Quad Shaping Drill 6mm(D) X 11.5mm(L)	QSD611	•	•

Description		Item #		
11	Quad Shaping Drill 6mm(D) X 13mm(L)	QSD613	•	•
11	Quad Shaping Drill 6mm(D) X 15mm(L)	QSD615	•	•
12	Narrow Posterior Hexed Driver - 17mm(L)	PHD00N	•	•
13	Posterior Large Hexed Driver - 17mm(L)	PHD02N	•	•
14	Narrow Right Angle Large Hexed Driver Tip - 24mm(L)	RASH3N	•	•
15	Certain Universal Ratchet Ext. Driver 4.1, 5, 6mm(D) - Short	IRE100U	•	•
16	Certain Universal Ratchet Ext. Driver 4.1, 5, 6mm(D) - Long	IRE200U	•	•
17	Direction Indicator 2mm(D) X 10mm(L) 3 Ea.	DI100	•	•
18	ACT Drill Extension	ACTDE	•	•
19	Handpiece Connector	MDR10	•	•
20	Implant Mount 3.4mm(D) X 15mm(L)	MMC15	•	•
21	Certain Implant Driver Tip 4, 5, 6mm(D) - Long	IIPDTL	•	•
22	Certain Implant Driver Tip 4, 5, 6mm(D) - Short	IIPDTS	•	•
23	Certain Implant Driver Tip 3.4mm(D) - Short	IMPDTS	•	•
24	Certain Implant Driver Tip 3.4mm(D) - Long	IMPDTL	•	•
25	Tapered Implant Dense Bone Tap 3.25mm(D) X 8.5mm(L)	NTAP3285	•	•
25	Tapered Implant Dense Bone Tap 3.25mm(D) X 10mm(L)	NTAP3210	•	•
25	Tapered Implant Dense Bone Tap 3.25mm(D) X 11.5mm(L)	NTAP3211	•	•
25	Tapered Implant Dense Bone Tap 3.25mm(D) X 13mm(L)	NTAP3213	•	•
25	Tapered Implant Dense Bone Tap 3.25mm(D) X 15mm(L)	NTAP3215	•	•
26	Tapered Implant Dense Bone Tap 4mm(D) X 8.5mm(L)	NTAP485	•	•
26	Tapered Implant Dense Bone Tap 4mm(D) X 10mm(L)	NTAP410	•	•
26	Tapered Implant Dense Bone Tap 4mm(D) X 11.5mm(L)	NTAP411	•	•
26	Tapered Implant Dense Bone Tap 4mm(D) X 13mm(L)	NTAP413	•	•
26	Tapered Implant Dense Bone Tap 4mm(D) X 15mm(L)	NTAP415	•	•
27	Tapered Implant Dense Bone Tap 5mm(D) X 8.5mm(L)	NTAP585	•	•
27	Tapered Implant Dense Bone Tap 5mm(D) X 10mm(L)	NTAP510	•	•
27	Tapered Implant Dense Bone Tap 5mm(D) X 11.5mm(L)	NTAP511	•	•
27	Tapered Implant Dense Bone Tap 5mm(D) X 13mm(L)	NTAP513	•	•
27	Tapered Implant Dense Bone Tap 5mm(D) X 15mm(L)	NTAP515	•	•
28	Tapered Implant Dense Bone Tap 6mm(D) X 8.5mm(L)	NTAP685	•	•
28	Tapered Implant Dense Bone Tap 6mm(D) X 10mm(L)	NTAP610	•	•
28	Tapered Implant Dense Bone Tap 6mm(D) X 11.5mm(L)	NTAP611	•	•
28	Tapered Implant Dense Bone Tap 6mm(D) X 13mm(L)	NTAP613	•	•
28	Tapered Implant Dense Bone Tap 6mm(D) X 15mm(L)	NTAP615	•	•
	Open End Wrench*	CW100	•	•
	Certain Driver Replacement O-Rings (10-Pack)*	IRORDR	•	•
	Tapered Implant Dense Bone Tap Kit Tray*	NTAPKT	•	•
	Ratchet Extension Short 6mm(L)*	RE100	•	•
	Ratchet Extension Long 15mm(L)*	RE200	•	•
	Ratchet Wrench*	WR150	•	•
	Base, Molded Insert & Smoked Lid*	QNTSKT	•	•
	Implant Depth Gauge 20mm(L)*	DPO20	•	•

* Items not pictured in this brochure.

Looking For Simplicity In Site Preparation? Choose The QNTSK Surgical Kit Today!

VISIT:
www.biomet3i.com

For More Information, Please Contact Your Local BIOMET 3i Sales Representative.



Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Toll Free: 1-800-342-5454
Phone: +1-561-776-6700
Fax: +1-561-776-1272
www.biomet3i.com



BIOMET 3i Dental Ibérica S.L.
EMEA Headquarters
WTC Almeda Park, Ed. 4, Planta 2
C/Tirso de Molina, 40
08940, Cornellà de Llobregat (Barcelona) Spain
Phone: +34-93-470-55-00
Fax: +34-93-371-78-49
www.biomet3i.com

ACT, Certain and PREVAIL are registered trademarks of BIOMET 3i LLC. Providing Solutions – One Patient At A Time is a trademark of BIOMET 3i LLC. ©2016 BIOMET 3i LLC.

All trademarks herein are the property of BIOMET 3i LLC unless otherwise indicated. This material is intended for clinicians only and is NOT intended for patient distribution. This material is not to be redistributed, duplicated, or disclosed without the express written consent of BIOMET 3i. For additional product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see the product package insert and the BIOMET 3i Website.



INST1096
REV B 02/16