


Want Simplicity
In Site Preparation?
Try The
Parallel Walled
Surgical Kit!



DIMENSIONS OF SUCCESS

The Paralled Walled Surgical Kit

Site Preparation Flowchart
For Drilling Simplicity

Color Coding For Ease Of Use



Convenient Pre-staging
Of Implants

Designed For All **Biomet 3i** Parallel
Walled Implants For Versatility


Providing Clinicians One Solution At A Time With The PSKT Surgical Kit

- Provides For Appropriate Instrument Selection And Simple Replacement Of Instrumentation
- Designated Grommets For ACT® And Countersink Drills
- Versatility Of One Surgical Kit For Multiple Implants
- Implant Diameter Illustrated From Wide To Small

BIOMET 3i™
PROVIDING SOLUTIONS - ONE PATIENT AT A TIME™

*For information on the drilling guidelines for parallel walled implants, please refer to the Surgical Manual (CATSM).

Simple Steps For The Selection



30
PHD00N



31
PHD02N



32
RASH3N



33
MTAP2



34
TAP413



35
TAP53S





36
TAP63S



37
DI100

And Replacement Of Instrumentation





OSSEOTITE[®] Dense Bone Tap Ordering Information



Length	5.0mm(D)	6.0mm(D)
6.5/7.0(L) - 8.5mm(L)	XTAP58S	XTAP68S
6.5/7.0(L) - 13.0mm(L)	XTAP53S	XTAP63S
6.5/7.0(L) - 18.0mm(L)	XTAP518S	XTAP618S

Note: Tapping in dense bone is required. Exceeding an insertion torque of more than 90Ncm may deform or strip the driver tip or the implant's internal hex and may possibly delay the surgical procedure.

PSKT Surgical Tray Contents

	Description	Item Number	PSKT40	PSKT35	PSKT30	PSKT30U	PSKT20	PSKT10	PSKT01
1	Round Drill	RD100							
2	Drill Extension	ACTDE							
3	ACT® Twist Drill 2.0mm X 10.0mm	ACT2010							
3	ACT® Twist Drill 2.0mm X 15.0mm	ACT2015							
3	ACT® Twist Drill 2.0mm X 20.0mm	ACT2020							
4	Pilot Drill 2.0/3.0mm	PD100							
5	ACT® Twist Drill 2.75mm X 10.0mm	ACT2710							
5	ACT® Twist Drill 2.75mm X 15.0mm	ACT2715							
5	ACT® Twist Drill 2.75mm X 20.0mm	ACT2720							
6	ACT® Twist Drill 3.0mm X 10.0mm	ACT3010							
6	ACT® Twist Drill 3.0mm X 15.0mm	ACT3015							
6	ACT® Twist Drill 3.0mm X 20.0mm	ACT3020							
7	Inner Surgical Kit Tray	PTT300UI							
8	ACT® Twist Drill 3.25mm X 10.0mm	ACT3210							
8	ACT® Twist Drill 3.25mm X 15.0mm	ACT3215							
8	ACT® Twist Drill 3.25mm X 20.0mm	ACT3220							
9	Countersink Drill 4.0mm (Ext. Con.)	CD100							
10	Certain® Countersink Drill 4.0mm	ICD100							
11	Countersink Depth Indicator For ICD100	ICDI100							
12	Countersink Drill for OSSEOTITE XP® 4.0/5.0mm	CD4500							
13	Countersink Depth Indicator For CD4500	CDI4500							
14	Pilot/Countersink Drill 5.0mm	CD500							
15	ACT® Twist Drill 3.85mm X 10.0mm	ACT3810							
15	ACT® Twist Drill 3.85mm X 15.0mm	ACT3815							
15	ACT® Twist Drill 3.85mm X 20.0mm	ACT3820							
16	ACT® Twist Drill 4.25mm X 8.5mm	ACT4285							
16	ACT® Twist Drill 4.25mm X 13.0mm	ACT4213							
16	ACT® Twist Drill 4.25mm X 18.0mm	ACT4218							
17	Countersink Drill for OSSEOTITE XP® 5.0/6.0mm	CD5600							
18	Countersink Depth Indicator For CD5600	CDI5600							
19	Pilot/Countersink Drill 6.0mm	CD600							
20	ACT® Twist Drill 4.85mm X 8.5mm	ACT4885							
20	ACT® Twist Drill 4.85mm X 13.0mm	ACT4813							
20	ACT® Twist Drill 4.85mm X 18.0mm	ACT4818							
21	ACT® Twist Drill 5.25mm X 8.5mm	ACT5285							
21	ACT® Twist Drill 5.25mm X 13.0mm	ACT5213							
21	ACT® Twist Drill 5.25mm X 18.0mm	ACT5218							
22	Certain® Ratchet Ext. Driver 3.4mm(D) - Short	IMRE100							
23	Certain® Ratchet Ext. Driver 3.4mm(D) - Long	IMRE200							
24	Certain® Implant Driver Tip 3.4mm(D) - Short	IMPDTS							
25	Certain® Implant Driver Tip 3.4mm(D) - Long	IMPDTL							
26	Certain® Ratchet Extension Driver 6.0mm(L)	IRE100							
27	Certain® Ratchet Extension Driver 15.0mm(L)	IRE200							
28	Certain® Implant Driver Tip - Short	IIPDTS							
29	Certain® Implant Driver Tip - Long	IIPDTL							
30	Narrow Small Hex Hand Driver	PHD00N							
31	Narrow Large Hex Hand Driver	PHD02N							
32	Narrow Large Hex Driver Tip	RASH3N							
33	Tap 3.25mm(D) X 18.0mm(L)	MTAP2							
34	Tap 4.0mm X 13.0mm	TAP413							
35	Tap 5.0mm X 13.0mm	TAP53S							
36	Tap 6.0mm X 13.0mm	TAP63S							
37	Direction Indicator 2.0mm X 10.0mm (3 each)	DI100							
38	Mount 3.4mm(D) X 15.0mm(L)	MMC15							
39	Handpiece Connector	MDR10							
40	Ratchet Extension 6.0mm	RE100							
41	Ratchet Extension 15.0mm	RE200							
42	OSSEOTITE® 2 Tap 5.0mm X 13.0mm	XTAP53S							
43	OSSEOTITE® 2 Tap 6.0mm X 13.0mm	XTAP63S							
44	Plastic Surgical Tray*	PTT300I							
45	Open End Wrench*	CW100							
46	Ratchet Wrench*	WR150							
47	Implant Depth Gauge*	DP020							
48	Plastic Organizer Tray*	PTT100							
49	Driver Tip Replacement O-Rings - 10 pack*	IRORDR							

*Items not pictured in this brochure

**For More Information On How Simple Site Preparation
Can Be With The Right Surgical Kit, Please Contact
Your **BIOMET 3i**™ Representative Today Or Visit**

www.biomet3i.com

BIOMET 3i
4555 Riverside Drive
Palm Beach Gardens, FL 33410
1-800-342-5454
Outside The U.S.: +1-561-776-6700
Fax: +1-561-776-1272
www.biomet3i.com

EC REP BIOMET 3i
Europe, Middle East & Africa
WTC Almeda Park, Ed. 1, Planta 1ª
Pl. de la Pau, s/n
08940, Cornellà de Llobregat
(Barcelona) Spain
Phone: +34-93-470-55-00
Fax: +34-93-371-78-49

Enhanced
Knowledge

Comprehensive
Treatment

Business
Excellence



ART979

REV E 08/11

Not Available In All Markets. Please Consult Your Local BIOMET 3i Sales Representative For Availability



**Scan With Your
Smartphone!**

In order to scan QR codes,
your mobile device
must have a QR code read-
er installed.

BIOMET 3i™
PROVIDING SOLUTIONS – ONE PATIENT AT A TIME™



Join Us



Follow Us



Watch Us

ACT, Certain and OSSEOTITE XP are registered trademarks and Providing Solutions – One Patient At A Time and design are trademarks of BIOMET 3i LLC. BIOMET 3i and design are trademarks and BIOMET is a registered trademark of BIOMET, Inc. ©2011 BIOMET 3i LLC. All rights reserved.