

# I-SERIES

compatible to: Biomet 3i® / Osseotite® External Hex®



## PLATFORM

3,4 mm

4,1 mm

5,0 mm

## TITANIUM BASE

for individual milled  
Zirconium Abutment  
incl. screw  
Titanium Grade 5 ELI



I 800

I 810

I 820



## SCAN BODY 3D-GUIDE®

for intra + extra  
oral scanning



I 9.S3D3.400

I 9.S3D4.150

I 9.S3D4.150



## LAB ANALOG

Stainless Steel



I 50

I 51

I 52



## STRAIGHT ABUTMENT

Gingiva Height 2,5 mm  
incl. screw  
Titanium Grade 5 ELI



—

I 110

—



## IMPLANT PICK-UP

for open-tray with  
integrated screw  
Coming soon



I TR-00013.4

I TR-00024.1

—

## ABUTMENT SCREW

Recommended tightening  
torques 35 Ncm  
Titanium Grade 5 ELI



I 61

I 61

I 61



Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA

# I-SERIES

compatible to Biomet 3i® /  
Osseotite® External Hex®



4,1 mm



I 10.Mi4.100

★

I 10.Mi4.100.05

★



I 5.La4.100

★



W 12.200.000

★



I 9.S3D4.150

★



I 6.Sc4.100

★



I 61

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA 

# I-SERIES DIM ANALOG

compatible to: Biomet 3i® / Osseotite® External Hex®



**PLATFORM**

3,4 mm

4,1 mm

5,0 mm

**DIM ANALOG**



I 5.DIM.340

I 5.DIM.410

I 5.DIM.500



**SLEEVE**

for plaster model  
Stainless Steel



DIM 5.XSL.A10

DIM 5.XSL.A20

DIM 5.XSL.A30



Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.  
Products indicated with ★ are available in the USA 