

L-SERIES

compatible to: Straumann / Bone Level®



PLATFORM

3,3 mm NC

4,1 / 4,8 mm RC

TITANIUM BASE

for individual milled
Zirconium Abutment
incl. screw
Titanium Grade 5 ELI



L 800



L 810

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



L 9.S3D3.300



L 9.S3D4.148

★

LAB ANALOG

Stainless Steel



L 50



L 51

★

IMPLANT PICK-UP

for open-tray with
integrated screw



L TR-NC013.3



L TR-RC024.1

ABUTMENT SCREW

Recommended tightening
torques 35 Ncm
Titanium Grade 5 ELI



L 60



L 60

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

L-SERIES

compatible to: Straumann / Bone Level®



PLATFORM

3,3 mm NC

4,1 / 4,8 mm RC

STRAIGHT ABUTMENT

Gingiva Height 1,5 mm

incl. screw

Titanium Grade 5 ELI



L 100



L 110

STRAIGHT ABUTMENT

Gingiva Height 3,0 mm

incl. screw

Titanium Grade 5 ELI



L 100-3



L 110-3

ABUTMENT SCREW

Recommended tightening

torques 35 Ncm

Titanium Grade 5 ELI



L 60



L 60



Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

L-SERIES

compatible to: Straumann / Bone Level®



PLATFORM

3,3 mm NC

4,1 / 4,8 mm RC

ANGLED ABUTMENT 18°

angled over surface
Gingiva Height 1,5 mm
incl. screw
Titanium Grade 5 ELI



L 200-1



L 210-1

ANGLED ABUTMENT 18°

angled over edge
Gingiva Height 1,5 mm
incl. screw
Titanium Grade 5 ELI



L 200-2



L 210-2

ANGLED ABUTMENT 18°

angled over surface
Gingiva Height 3,0 mm
incl. screw
Titanium Grade 5 ELI



L 200-1-3



L 210-1-3

ANGLED ABUTMENT 18°

angled over edge
Gingiva Height 3,0 mm
incl. screw
Titanium Grade 5 ELI



L 200-2-3



L 210-2-3

ABUTMENT SCREW

Recommended tightening
torques 35 Ncm
Titanium Grade 5 ELI



L 60



L 60



Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

L-SERIES

compatible to: Straumann / Bone Level®



PLATFORM

3,3 mm NC

4,1 / 4,8 mm RC

2-CONNECT®-ABUTMENT SET

incl. 2-CONnect®-Coping
and Coping-screw
Titanium Grade 5 ELI



Gingiva Height ... 0,5 mm
... 2,5 mm

L 800 S

L 810 S



L-2CON 8033002S

L-2CON 8041482S



2-CONNECT®-ABUTMENT

Recommended tightening
torques 35 Ncm

Titanium Grade 5 ELI

M (Male Part)



Gingiva Height ... 0,5 mm
... 2,5 mm

L 800 M

L 810 M



L-2CON 8033002M

L-2CON 8041482M



2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U



COPING-SCREW

Recommended tightening
torques 30 Ncm

Titanium Grade 5 ELI



KS 61

KS 61



Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 

L-SERIES

compatible to: Straumann / Bone Level®



PLATFORM

3,3 mm NC

4,1 / 4,8 mm RC

2-CONNECT®-ABUTMENT SET

incl. 2-CONnect®-Coping,
Coping-screw & Male Part Screw
angled 20°

Gingiva Height 2,5 mm

Titanium Grade 5 ELI

Coming soon



over surface

L-2CON 80330A3S

L-2CON 84148A3S

over edge

L-2CON 80330C3S

L-2CON 84148C3S

2-CONNECT®-ABUTMENT

angled 20°

Gingiva Height 2,5 mm

Titanium Grade 5 ELI

Coming soon

M (Male Part)



over surface

L-2CON 80330A3M

L-2CON 84148A3M

over edge

L-2CON 80330C3M

L-2CON 84148C3M

2-CONNECT®-COPING

Titanium Grade 5 ELI



F (Female Part)

2CON F-U

2CON F-U



COPING-SCREW

Recommended tightening

torques: 15 Ncm

Titanium Grade 5 ELI

Coming soon



KS 62

KS 62

MALE PART SCREW

for angled 20° 2-CONnect®

Abutment | Recommended

tightening torques: 20 Ncm

Titanium Grade 5 ELI

Coming soon



L 6.Ms3.341

L 6.Ms3.341

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

L-SERIES

compatible to: Straumann / Bone Level®



PLATFORM

3,3 mm NC

4,1 / 4,8 mm RC

2-CONNECT®- IMPLANT PICK-UP

integrated screw



● Straight KS 61 = $\varnothing$ 1,8 mm

2CON TR-00021.8

2CON TR-00021.8

● Angled KS 62 = $\varnothing$ 1,6 mm

2CON TR-00031.6

2CON TR-00031.6

2-CONNECT®- SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



● Straight KS 61 = $\varnothing$ 1,8 mm

2CON 9.S3D1.800

2CON 9.S3D1.800

★

● Angled KS 62 = $\varnothing$ 1,6 mm

2CON 9.S3D1.600

2CON 9.S3D1.600

★

2-CONNECT®- LAB ANALOG

Stainless Steel



● Straight KS 61 = $\varnothing$ 1,8 mm

2CON 5.La1.800

2CON 5.La1.800

★

● Angled KS 62 = $\varnothing$ 1,6 mm

2CON 5.La1.600

2CON 5.La1.600

★

2-CONNECT®- DIM ANALOG



● Straight KS 61 = $\varnothing$ 1,8 mm

2CON 5.DIM.180

2CON 5.DIM.180

● Angled KS 62 = $\varnothing$ 1,6 mm

2CON 5.DIM.160

2CON 5.DIM.160

★

SLEEVE

plaster model
Stainless Steel



DIM 5.XSL.A30

DIM 5.XSL.A30

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 

L-SERIES

compatible to: Straumann / Bone Level®



PLATFORM

3,3 mm NC

4,1 / 4,8 mm RC

NT-PREFORM®

incl. coloured milling screw
Titanium Grade 5 ELI



Gingiva Height

1,3 mm	L 10.Mi3.300	L 10.Mi4.148
0,5 mm	L 10.Mi3.300.05	L 10.Mi4.148.05

NT-MILL LAB ANALOG

for nt-Preform®
Stainless Steel



L 5.La3.300

L 5.La4.148

★

ABUTMENT GRIP

for nt-Holder fixture
Stainless Steel



W 12.100.000

W 12.200.000

★

SCAN BODY 3D-GUIDE®

for intra + extra
oral scanning



L 9.S3D3.300

L 9.S3D4.148

★

COLOURED MILLING SCREW

for nt-Preform®
Recommended tightening
torques 10 Ncm



L 6.Sc3.341

L 6.Sc3.341

★

ABUTMENT SCREW

Recommended
tightening torques
35 Ncm
Titanium Grade 5 ELI



L 60

L 60

★

Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA

L-SERIES DIM ANALOG

compatible to: Straumann / Bone Level®



PLATFORM

3,3 mm NC

4,1 / 4,8 mm RC

DIM ANALOG



L 5.DIM.330

L 5.DIM.414



SLEEVE

for plaster model
Stainless Steel



DIM 5.XSL.A10

DIM 5.XSL.A20



Notice: Products indicated with ® and ™ are registered brand names of the respective manufactures.
Products indicated with ★ are available in the USA 